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The screenshot shows a software interface for protein detection. The main title is '蛋白质检测' (Protein Detection). There are several sections:

- 人员管理 (Personnel Management):** Fields for '检测人员' (Inspector), '被检单位' (Sample Unit), '检测单位' (Testing Unit), and '联系电话' (Contact Number).
- 检测结果 (Test Results):** Displays '滴定体积' (Titration Volume) as 8.621 and '含氮量' (Nitrogen Content) as 6.1267%. Below it, '蛋白质含量' (Protein Content) is shown as 38.2920%.
- 项目管理 (Project Management):** Lists parameters for the current test: '项目名称' (Project Name) is 大豆蛋白, '标准酸浓度(mol/L)' (Standard Acid Concentration) is 0.1059, '样品质量(g)' (Sample Mass) is 0.2046, '空白体积(mL)' (Blank Volume) is 0.166, '加酸体积(mL)' (Add Acid Volume) is 50, '加稀释液体积(mL)' (Add Diluent Volume) is 20, '加碱体积(mL)' (Add Alkali Volume) is 60, '蒸馏时间(min)' (Distillation Time) is 6, '蛋白质系数' (Protein Coefficient) is 6.25, and '指示剂' (Indicator) is 甲基红:溴甲酚绿=1:5.
- 操作按钮 (Operation Buttons):** Includes '自动检测' (Automatic Detection), '空白检测' (Blank Detection), '样品检测' (Sample Detection), and '手动检测' (Manual Detection).
- 底部功能栏 (Bottom Function Bar):** Contains buttons for '项目类型' (Project Type), '保存项目' (Save Project), '设置系数' (Set Coefficient), '保存人员' (Save Personnel), '打印' (Print), '上传' (Upload), and '仪器调试' (Instrument Adjustment).

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